

Basic User Guide

Guide to Basic OnePacs Functionality

Guide to the basic steps involved in study lookup and report review



- [Guide to Basic OnePacs Functionality](#)
- [Purpose and Scope](#)
- [Browser Support](#)
- [Logging in](#)
 - [Resetting your password](#)
 - [2-factor authentication](#)
- [Worklist](#)
 - [Worklist review](#)
 - [Case data displayed on worklist](#)
 - [Worklist Navigation](#)
 - [Viewing and printing reports](#)
- [Additional features of OnePacs](#)

Purpose and Scope

The purpose of this brief guide is to familiarize the new user of OnePacs with the steps involved in searching for studies, reviewing study status, using the simple JPEG based image viewer, and reviewing reports.

There are many features available in OnePacs which are not covered in this Quick-start guide. For more detailed information on using and configuring OnePacs, please consult your local system administrator or your designated OnePacs customer support representative.

Browser Support

OnePacs supports the following browsers:

- Chrome 52+ (recommended)
- Firefox 49+
- Internet Explorer 11 with latest updates (some features may be unavailable)

Logging in

First, connect to the website specified in your welcome email.

This will redirect to a secure, encrypted login page. You will then be presented with the login screen shown in the figure below.

Type your user name in the text field marked "User", and your password in the text field marked "Password". When this is done, click "Login".



Login @onepacs.com

User Id:
Password:



[Forgot your password?](#)

Resetting your password

If you have forgotten your password click the "Forgot your password" link and complete the form on the next page. A reset password link will be emailed to you or to your OnePacs administrator.

2-factor authentication

If your user account is configured to require 2-factor authentication, you will be prompted to enter a 6-digit security code generated from a registered mobile device. Enter the code and click "Login" to continue.

Login @onepacs.com

Security code:



[Forgot your password?](#)



2-factor Authentication

2-factor authentication adds an additional layer of security to your OnePacs account and is recommended for all users. To add 2-factor authentication to your account, navigate to Options -> User Preferences and click on the "2-factor Authentication" button.

Worklist

The "worklist" functionality of OnePacs refers to the standard means of listing cases which are in the OnePacs system according to the user's desired criteria. Users have the option of generating one or more customized worklist views, which can be accessed through tabs at the top of the OnePacs worklist view page. The two main types of worklist views are "search views", which dynamically search the system for cases that meet certain criteria, and "folders", which are static lists of cases that are placed into folders manually by the user (or, which can be auto-added to the folder by OnePacs according to auto-store criteria). When a case is added to a folder, it remains in the folder until manually removed. This type of worklist view is appropriate for management of teaching files, or for storing cases which need to be acted on later by the radiologist (such as cases requiring followup or phone calls).

In contrast, search views are continuously updated according to the user's chosen criteria, and cases which no longer meet the criteria for the search view will be removed from the worklist. Search views are useful for workflow management, in that cases which no longer meet appropriate criteria (such as, "unread cases from the last 3 days") will be removed from the worklist when their status changes (such as, when the case status changes from "unread" to "read").

A more detailed overview of the various components of the worklist view may be found [here](#).

Worklist review

After you log in, the system will display a list of studies and information about those studies. By default, the system will list all studies which have been sent to the OnePacs system within the last 30 days. This time range is indicated in the upper right hand corner, in the box which contains the text "Last 30 days". (See figure below.) Criteria that results in studies being filtered will be highlighted in yellow. Click the green arrows on the search panel to reset the search criteria.

Toolbar	Arrive Time	TT...	Study Date	Name	Birth Date	MRRN	Body Part/Description	History	Arrived	Age	Mod	Facility	Assigned to	Physician	Img
	09/11 10:46 CDT	-89d	09/11 10:45 CDT	Sample Case		SAMPLECASE	CHEST	pt history (HX)	89d	DX	Demo Facility A	Falk, Justin	Dr Demo	2	
	09/11 10:46 CDT	-89d	09/11 10:45 CDT	Sample Case		SAMPLECASE	RIGHT WRIST	test (HX)	89d	MR	Demo Facility A	Admin, Demo	Admin, Demo	59	
	09/11 10:46 CDT	-89d	09/11 10:45 CDT	Sample Case		SAMPLECASE	BRAIN	test (HX)	89d	MR...	Demo Facility A	Admin, Demo	Admin, Demo	141	
	10/13 11:51 CDT	-57d	10/13 11:49 CDT	Sample Case		SAMPLECASE	CHEST		57d	CT	Demo Facility C	Admin, Demo	Admin, Demo	111	
	10/27 16:19 CDT		10/13 11:49 CDT	Sample Case		SAMPLECASE	ESOPHAGUS		43d	RF	Demo Facility B	Admin, Demo	Admin, Demo	42	
	09/11 10:46 CDT		09/11 10:45 CDT	Sample Case		SAMPLECASE	RIGHT WRIST	test (HX)	89d	MR	Demo Facility C	Falk, Justin	Dr. John Doe	59	
	10/14 16:39 CDT		10/13 11:49 CDT	Sample Case		SAMPLECASE	BRAIN		57d	MR	Demo Facility C	Admin, Demo	Admin, Demo	139	
	09/11 10:46 CDT		09/11 10:45 CDT	Sample Case		SAMPLECASE	CHEST	mva-Lt. side chest & abd. pain. (HX)	89d	CT	Demo Facility A	Admin, Demo	Admin, Demo	111	
	10/13 11:51 CDT		10/13 11:49 CDT	Sample Case		SAMPLECASE	CHEST		57d	DX...	Demo Facility C	Admin, Demo	Jane Smith, MD	3	
	09/05 15:29 CDT		09/05 15:28 CDT	Sample Case		SAMPLECASE	PELVIS		95d	CR	Demo Facility A	Radiologist, Demo		2	

Case data displayed on worklist

Each column of the worklist shows information about one study. The information displayed includes the patient's name, the date the study was performed, the date the study arrived on the OnePacs system, as well as the current status of the study.

Study status icons

Studies will have one of the following statuses.

- In transmission to OnePacs: the status column will show an hourglass icon:
- Arrived at OnePacs, but awaiting confirmation of the study by the transmitting facility: A blank report page with a question mark will display:
- Arrived at OnePacs and ready for a radiologist to read: a blank report page icon will display:
- Disregarded – if a study was received by OnePacs but will not receive a report (examples would include studies sent only for comparison purposes), the blank report icon will display, with a red "minus sign" symbol overlaid on the report:
- If a report is currently being written for a study, a lock icon will display:
- After a study has been read and the report for the study is available, the report paper icon will show lines of text and a green check mark:

For studies which have reports, a second icon will appear on the worklist which is used to view the radiology report. This icon consists of a pair of eyeglasses overlaid on a piece of paper: This icon will appear immediately on the right of the "report ready icon", as follows:

Thus, on the sample set of cases seen in the figure below,

The 1st case is currently being sent to OnePacs.

The 2nd case has arrived at OnePacs, but is awaiting confirmation of the study order and details.

The 3rd study has been confirmed for reading and is locked for reporting.

The 4th case has been disregarded.

The 5th study is ready for a radiologist to read, but have not been read or assigned yet.

The 6th study was sent only for comparison and does not require a report.

The 7th study has a report, which is available to read by clicking the display report icon.

Toolbar	Send Time	Name
	06/14 9:17	SAMPLE CASE,
	06/14 9:09	SAMPLE CASE,


Study statuses are shown in the left hand column.

Worklist Navigation

Paging Controls

The worklist may display more cases than fit on one screen. In this case, the list of cases will be split across more than one page. At the bottom left hand corner of the worklist view, there are navigation arrows that allow the user to switch between “pages” of a multiple-page list of cases. (See figure below.)

The paging controls will generally indicate how many pages of cases there are, and which page is currently being displayed. If there are more than 1000 results some paging controls may not be available and the total results will only show that there is more than 1000 results. The arrow pointing to the right can be clicked on to switch to the next page of cases, and the arrow pointing to the left can be clicked to switch to the previous page of cases. Text will pop-up indicating which arrow switches to the previous page, and which to the next page of cases. If there are more than 1000 results and you want to see the total number of results for the current worklist filter you can hover over the estimated results to query for the total number of results.





Paging controls. At the lower left hand corner of the screen, controls exist to switch between multiple pages of cases when there are more cases than fit on a single screen.

Date range

The time period can be changed to different criteria. For example, for night call one may wish to see only cases from the previous 6 hours. Or, to look up results for an old case, one may wish to see cases from the last 30 days. By selecting the appropriate date range, one can change the set of cases which display on the screen to facilitate finding the cases one needs. (See figure below.). Various date search criteria are also available.

Date ranges




Study Status checkboxes

In the upper right hand corner of the screen, there is a list of study statuses. The checkboxes indicate whether you wish for studies that have that status to appear on the study list. When you first log in, the system will display cases of any status, therefore all the checkboxes are checked. If you are only interested in seeing cases that have already been read already, you can de-select all of the checkboxes except the "Read" button. (See figure below.)



Study status checkboxes, with all boxes checked. The worklist will show all cases, of any status.




Study status checkboxes, with only the "Read" box checked. The worklist will show only cases that have already been read, and have reports available in any confirmation state.

Searching for cases

After setting the date range, and selecting the checkboxes for the types of statuses for which you want studies to display on the worklist, you can next narrow down the cases that display on the list by entering in search criteria. In the upper portion of the screen, there are several search boxes. In Figure 7, the "Patient Name" and "Study Description" search boxes are circled.

To limit studies displayed on the worklist by search criteria, simply enter your desired search text in the text field under the listed search criteria. For example, one may search for a specific patient with a certain study description.



Viewing and printing reports

When a report has the status of "complete", there will be a view report icon on the worklist toolbar which appears displays the following icon in the toolbar:



When this icon is clicked, the report for the study will display in a new window (See figure below.) In the upper left corner of the window, there are icons that can be used to export the report to PDF, RTF (Microsoft Word), Print, or to view the images of the study in various viewers.

Example report display.



Report Submission Date: Jan 31, 2017 12:41:30 AM CST

Patient		Study	
Name:	SAMPLE CASE	Date:	Nov 26, 2016 5:49:36 AM CST
MRN:	SAMPLECASE	Modality Type:	RF
Gender:		Description:	ESOPHAGUS
DOB:		Body Part:	
Institution:	Main Street Hospital	DICOM Inst.:	SAMPLE CASE

History Rule out esophagitis

Comparisons:
None

Study: ESOPHAGUS

Contrast:

FINDINGS:

Lungs: Normal
Views:

Heart and Great Vessels: Normal

Mediastinum and Hila: Normal

Upper Abdomen: Normal

Additional abnormalities: None

Additional features of OnePacs

The following features of OnePacs may also be of interest:

Confirming studies

Depending on how OnePacs has been configured for your facility, it may be necessary for studies transmitted to OnePacs to be "confirmed" before they are ready to be interpreted for reporting. "Confirmation" of studies serves all the usual purposes of paper study requisitions - they confirm that interpretation has been requested for an examination, provide basic information regarding the patient and the requesting physician, describe where the report should be sent, and provide study history.

A full description of how to confirm studies is available [here](#).

Study Notes

Users with the View/edit study notes permission may add study related notes or flag a study. A full description is available [here](#)

Preferences

Users who have been given the "View user preferences" permission may access and modify their user settings, including user password.

To access the user preferences page, select the "User preferences" item from the "Options" menu:

For more information on User Preferences see [User Preferences](#).

Uploading images and CDs

CDs and other DICOM images may be uploaded directly from a facility. For more information see [Direct upload using the web interface](#)

Using the lightweight JPEG-based image viewer

A lightweight JPEG-based study viewer is available in OnePacs. Please [click here](#) for details.

Using the Quality Assurance (QA) case tracking system

OnePacs features a Quality Assurance (QA) case tracking system that enables authorized users to enter discrepancy reports for interpreted cases. This may be used by client radiology groups using the services of teleradiology groups, or any other individuals who may need to report discrepancies regarding examination interpretations.

The QA features facilitate keeping track of specific cases for which discrepancies need to be reported, allowing the interpreting radiologist and/or radiology group to review the QA issue and issue a response, and keeping track of statistics regarding QA issues.

Usage of the QA system is described in detail [here](#).